


Dental Led Curing Light Using Manual

1. Function button instruction



Intensity key: adjust luminance

Power switch: ON/OFF,lighting ON/OFF

Mode key: F0,F1,F2,F3. four working modes

Time setting key: Corresponding mode time selection, adjustment

Display area: time indicator

2.Function Instruction :

a) Press the power switch,the machine be standby state.

b) Mode key: F0,F1,F2,F3 four working modes,

F0 mode: 1 second mode

During the working time 1-5S, you can choose 1s 3s or 5s,luminance output 100% (can not be adjustment)

F1 mode: strong mode

During the working time 5-40S , the machine is lighting with high power.

you can press the side key “+ -” to adjust luminance,40%-100% can be adjustment,the luminance will be +1%/-1% per operation.

F2 mode: gradually strong mode

During the working time 5-40S , the machine light power is gradually strong in the first 10 seconds,and then remaining time is

lighting with constant high power . You can press the side key “+ -” to adjust luminance,40%-100% can be adjustment,the luminance will be +1%/-1% per operation.

F3 mode: flashing mode

During the working time 5-40s,the machine works in falshing mode with high power.you can press the side key “+ -” to adjust luminance,40%-100% can be adjustment,the luminance will be +1%/-1% per operation.

c) Press the power switch to stop working,the machine be standby state.

3:Technical data

Input voltage-----AC 110V—240V50—60Hz

Output voltage-----DC5V 1000mA

Charging voltage-----DC4.2V

Battery capacity-----2200mAh

wavelength-----455—480nm

Intensity -----1200—2600mw/cm²

4. Symbol Mark Instructions



part

attention,consult accompanying document



type B applied



向上

shipping package should be straight up during transportation



怕雨

keep dry



怕晒

used indoor only



class II equipment



Deal with the product according to the WEEE(2009/96/EC)



易碎物品 handle with care



Recovery




Storage temperature at -20°C~+55°C



Storage relative

humidity less than 80%



Storage atmospheric pressure at 70kpa~106kpa