

AIXIN MEDICAL

USER MANUAL

AX-J5

Polishing Lathe

with Suction System


Description:

The AX-J5 Polishing Lathe with Suction Box is used for polishing castings in dental laboratories, e.g., alloys (Co-Cr, Ni-Cr) stainless steel, flexible dentures, bridges and porcelain teeth. It features low noise, high rotational speed, built-in dust collector, safety devices, and easy operations. It has pass the CCMS safety certification.

Specifications:

Power	220V/110V 50/60Hz
Power Consumption	1,120W
Rotational Speed	3,000 RPM
Dimensions	50x20x25 (cm)

Diagram:



1. Suction Power
2. Light
3. Pan
4. Lathe

5. Protective shield
6. Light switch
7. Suction switch
8. Lathe switch

Instructions:

1. Unpack the package, check if all the accessories are included in the box, read the manual carefully before operating.
2. Check if the screws are loose and tighten them.
3. Rotate the lathe axis before operating, if there is friction please check and solve the problems before operating.
4. Install the protective shields, and adjust to the right position.
5. Fasten the base to the bench top if necessary.
6. Turn on the power to start operating the lathe.

Cautions:

1. Use the protective shields during polishing.
2. The lathe should be stopped for cooling after using for 30 minutes in order to protect the lathe.
3. Make sure the polishing wheels are tightly fixed to the spindle of the lathe.

Accessories:

User Manual	1
Buffing wheel	2 PCS