





Instructions for the Workstation to Open the Box and Install

- After removing the container, remove the relevant parts and check the packing list.
- I. Insert the flat tube ①(the position and direction of the screw hole is the same as the side plate) into the flat hole in the middle of the side plate, make sure the full table height is 80 or 85 cm, then tighten with the M8 hexagonal screw.
- II. Fasten the cast aluminum table foot ②with M8 inner hexagonal screw and the flat tube ①, measure the whole table height correctly, and then erect the whole table.
- III. Connect the lead wire of the bracket ③column to the lead wire on the desktop, then insert the rack post into the aluminum seat ④ of the desktop, tighten the sides with fixed screws respectively, and place the table in the reserved position. Connect the power line air source pipe into the electrical source assembly ⑤ (at the lower left of the table sandwich) and adjust the whole table by using the four corner adjusting screw ⑥ of the table foot.
- IV. Insert the working lamp ⑦ into the socket of the side of the post with wires in the desk hole⑧.
- V. Insert the elbow guard bracket ⑨ into the corresponding hole position ⑩ (remove the hole position panel to adjust the width of the elbow bracket for 40~48cm, factory set to 40cm)
- VI. Plug the vacuum assembly 11 into the flat hole 12, and the internal hose is connected. Press the vacuum spring plate 13 and insert the glass dust screen.
- VII. The vacuum cleaner power plug 14 is inserted into the power outlet 15 of the body to connect the body with the flat suction hose 16, and the flat suction other side interface is connected with the suction port hose.
- VIII. Pull out the trachea and connect it with the air gun 17 (if the trachea is too long, a section can be reduced appropriately, so that the air gun does not fall out after reduction)
- IX. After installation, turn on the power switch to confirm normal and start work.



- | | |
|------------------------------|------------------------------|
| 1-flat tube | 2-cast aluminum table foot |
| 3-bracket | 4-aluminum seat |
| 5-electrical source assembly | 6-adjustment screw |
| 7-working lamp | 8-desk hole |
| 9-elbow guard bracket | 10-hole position |
| 11-vacuum assembly | 12-flat hole |
| 13-vacuum spring plate | 14-vacuum cleaner power plug |
| 15-power outlet of the body | 16-flat suction hose |
| 17- air gun | |



(The diagram is 1. 2 m single mechanic table. 1. 0 m single table, 1. 4 m single table, 1. 6 m double table and 1. 8 m double table installation are the same.)