

OPERATION MANUAL

LED Curing Light

(Model: TR-DCL05)



Content

- 1 Introduction
- 2 Structure
- 3 Operation Instruction
- 4 Parameters
- 5 Precautions
- 6 Packing List

6. Packing List

- 1. LED Curing Light :1pc
- 2. Lamp Guide :1pc
- 3. Battery :1pc
- 4. Charger Socket :1pc
- 5. Adapter :1pc
- 6. Operation Manual :1pc
- 7. Light Protector :1pc



5. Precautions

1. Please change a new lamp guide if it is broken.
2. Disinfection: put the lamp guide to 75% alcohol for a few minutes.
3. Remove adapter from the power supply after charge.



4. Do not touch the battery electrodes by any metal objects, otherwise the battery will be damaged by the short-circuit.



5. It must be charged if it buzzs 5 times, otherwise the battery will be damaged because of the low power.




6. When you switch on the light and it switch off automatically, it only can be used after 1 min due to the overheat protection system.

1. Introduction :

Thank you for choosing our LED Curing Light JAS-2001B, with industrial intelligent chip, digital display and multi-functions, safe operation and perfect self-protection system. It has pure blue visible light and the wavelength is 400-490nm. It can cure all the dental resin in the market nowadays.

2. Structure



3. Operation Instruction

A. Charge

Put the curing light into the charger socket, the yellow light is on.

When it is done, the green light is on. It takes about 4-6 hours.

Remove the adapter from the power supply after charging.

B. Operation Details

Press the **ON/OFF** button, it is on.

Press the **ON/OFF** button over 3 seconds, it is off.

When it is on, you can set the TIME (from 1 - 6 seconds), then press the **ON/OFF** button, it is working; press the **ON/OFF** button again, it stops.

It will turn off automatically without any operation in a few minutes.

If it buzzes 5 times, it must be charged before using.

4. Parameters

Gross Weight : 576g

Input Voltage : AC 100-240V, 50/60Hz

Charger Socket : DC5V

Battery : 3500mAh

Wavelength : 385nm - 515nm

Light Density : $\geq 2700 \text{mw/cm}^2$.