

Contents

1. General information	2
2. Safety information	3
3. Construction brief instruction.....	4
4. Specification.....	6
5. Installation and Settings.....	7
6. Operation direction.....	9
7. Maintenance notice.....	18
8. Clearing of faults.....	20
9. Waste disposal.....	21
10. Taboos and precautions.....	22
11. Electrical schematic diagram of dental unit.....	24
12. Water and air connection diagram.....	25

1.General notice

1.1 Esteemed users

IT is our pleasure that you bought this dental unit made , this user manual carefully in order to master this machine as soon as possible

1.2 Usage notice

Read this user manual carefully in order to familiarize yourself with this new device before using remand keep in mind of the warning and safety notices listed in this book plate this user manual at a certain location where you can take easily to ensure that you and other dentists can get it handly.

1.3 warranty and duty

This device must be maintained regularly to ensure its safety and reliability in operation.thus to safeguard dentistsand patients or other'health and safety.

Maintenance can only be done by us or the authorised our agents,if there are something damaged with components which matter with safe performance of the device,the faulted parts should be chaged into our components.

1.4 Technician and the working area

Operators should be trained by our company or our agents,this device should be used for dental diagnosis and treatment,it can not used at the area with risk of explosion.

1.5 Accessories

The accessories mentioned in this manual like high speed handpiece,low-speed handpiece,scaler,curing light,are not included in the accesories of dental unit.they are the extra accesories available for further purchase.

2.2.Safety notice

2.1 Dangerous class of symbols

Comply with the warning and safety notice, to avoid casualty or property loss. please be attentioned to the marks below:

Warning
A possibly dangerous situation that could result in serious bodily injury or death.
Notice
A possibly harmful situation which could lead to damage of the product or an object in its environment.

Helpful information
Application instruction and other important information.

2.2 Symbol explanation

2.2.1 the labels as belows:




look to check with the enclosed document



B sytle applied parts

2.2.2 Explanation to the symbol on the wood boxes



No heavy pressure Avoid raining Upward placement, Handle with care keep out of sun

2.3 Media quality

To safeguard patients, the water supplied for this device must meet to request of drinking and the air should meet the request of medical level, it must be dry and oil-free.

To ensure requirment in water, user should schedully check the water quality.

2.4 Patient Chair

Maximum load capacity of the base of patient chair is 200KG.(treating devices on this dental unit total 245KG) patients must stretch feet and hands on the cushion.

2.5 Ventilation Slots

Under no circumstances may the ventilation slots on the unit

To over heat. Do not spray liquid such as disinfectants in into the ventilation slots. This may lead to malfunctions.

3.Construction

3.1 Construction diagram



R7 CONSTRUCTION

A	Base	D	Foot controller	G	Spittoon	J	Balance arm	M	Air lock
B	Cushion	E	Headrest	H	Operating light	K	Air pressure meter	N	Instrument rack
C	Backrest	F	Assistant control	I	Light arm	L	Touch sense Instrument tray		

3.2 Description of construction

This dental unit consists of patient chair, connector of handpiece, dental lamp, 3-way Syringe, weak saliva ejector, strong saliva ejector, spiton, viewer, foot controller, air and water supply system.

Electricity, Water and air are supplied from the chair base. There are air filter, Pressure reducing valve, water filter and electricity filter for purifying the water, air and electricity supplying to this device.

Air supply to handpieces is switched on or off with the foot controller. The air and water system of handpieces, ejectors are controlled by the relevant valves.

This dental unit is equipped with bottles for purified water, with which is helpful for handpieces to prolong its lifespan. Its design is easy for exchanging between tap water and purified water.

3.3 Application scation scope

Bing Used only for dental diagnosis and treatment.

4.Specification

Power supply: 230V AC \pm 22V , 50HZ \pm 1HZ;

Power consumption: \leq 600W;

Rated current: 3.2A;

Fuse: 05 \times 20, 6.3A;

Outside fuse for power supply:230V AC 16A (delaying style)

Safety class:Class II

Protection level of anti electric shock: A style applied parts

Work ways: Intermittent working

Signals connection: Refer to the parts with extra instruction

Signals connection class: Common electrical appliance;the foot control is for IPX4 level

Transportation and storage conditions:

Temperature:-40 $^{\circ}$ C--+55 $^{\circ}$ C

Relative humidity: 10%--90%

Air pressure :50kPa--106kPa

Working condition:Working temperature: 5 $^{\circ}$ C--40 $^{\circ}$ C

Air relative humidity: \leq 80%

No condensation

Air discharge volume: 55L/min;

Air input pressure: 500kPa--800kPa;

Water input pressure: 200kPa--400kPa;

Total weight: 245kg \pm 5kg;

Patient seat lowest height: \leq 410mm;

Patient seat highest height \geq 760mm;

Backrest front rake:105 $^{\circ}$ \pm 5 $^{\circ}$;

Back rake of Backrest: 180 $^{\circ}$ \pm 5 $^{\circ}$;

Air pressure for high speed handpiece: 0--300kPa;

Air pressure for high low speed handpiece:0--250kPa;

5.Installation and setting

5.1 Installation conditions

To make sure dental unit work well the conditions of air, water, electricity for the device should be


- Air supply air pressure: 0.5 MPa--0.8MPa Folw volume >50L/min
- Water supply Water pressure: 0.2MPa--0.4MPa Current volume>10Lmin
- Power supply AC 220V ± 10% 50Hz ± 1Hz

5.2 Preparation for installation

Decide the suitable location of installation place according to the elements loke general situation, daylighting,convenience for use. Where should be clean,dry,well ventilated,cool,and the floor must be solid,firm and flat.what’s more,the installation room should be wide enough for the arm stretch. After decising the location, you can connect the wires and tubes from outside ends as below:

5.2.1Arrangement of connecting water in/out,electricity,compressure air supply

Get the things below ready before installation A)End of power wire should stretch over 80mm from ground. B)Adopt the pipe of water or air with width 22 mm C)Outlet pipe of water should have internal diameter 40mm, stretch over 50mm from ground



5.2.2 Steps of installation

A) Patient chair installation

Fit the chair base at the planned place. If the installation floor is not flat enough, user can fix the five pcs of enclosed leveling screws (M12*15) into the five screwholes (M12), adjust the height of five screws to make sure the dental unit stand firmly onto the ground, and then fix four pieces of expansive screws into the ground, to make sure the device works stably and quietly.

B) Mount the patient chair cushion

Stretch the backrest to the end first, fix two flat head screws into screwholes of the seat through the screw slots of the base frame, adjust the seat board to best location to make sure the bottom surface touch firmly with the metal frame. Finally fix the rest two screws firmly.

C) Outlet unit installation

Seal this device's outlet tubes with the inlet tubes from outside tightly, thus to avoid the water and air input or output leak outside. The outlet unit may set in the frame base or outside of the base, it depends on your choice. Open the cover you can find fix the outlet tubes.

D) Tubes connection

You must get rid of the water and air inside this dental unit, as well as the waste, note to get the tubes connected tightly with outside tubes.

E) Dental lamp installation

Pass the power wire through the supporting arm inside, and connect it with the wire from pillar bases of the dental unit, put the clasp end into the pillar and then fix the supporting pillar onto the dental unit, (mind not to damage the power wire when connecting wire)

Pass the power wire from the pillar into the cover ring of dental to fit together with the wire of dental lamp. Insert the connective pivot into the balance arm, and fix with screw.

F) Connect power supply

Put the plug into the power socket, which must have the wire connected to ground. The wire should have the copper wire inside less than 1mm, make sure this wire connected well with ground.

5.3 Modulation

Now the device is at state of modulation after different parts are fitted well.

Open the base box, check the air pressure meter of air filtering adjusting valve whether it reads 0.6mpa, if not, adjust the valve to keep the air pressure reading 0.6mpa. (Method: pull up the rotary knob at the top air pressure valve, turn clockwise for increasing pressure, turn anticlockwise for reducing pressure.)

6.Operation instruction

6.1 Master switch of water,electricity,and air

This switch is set at the back of base,Pulling out means power,water and air on,Pushing down means power,water and air off.

Notice

For the sake of protection,prolong lifespan of the device, as well as saving power,please get this mater switch off.

6.2 Dental lamp

Dental lamps are devided into : Led dental lamp

6.2.1 LED dental lamp

LED dental lamp have a function of continuously adjusting classless luminance. Please refer to"LED Dental Lamp User Manual". Power on or off can be controlled by foot controller or main control panel. Press the button down first to start. Press again to shut down the light. Every pressing will receive buzzer feed-back.

Warning

No shining eyes directly with this LED dental lamp

6.3 Connective end for high speed handpiece

Take the connective end of high speed handpiece off the instrument rack,connect it to high-speed handpiece.step on the side of foot controller,which mark"W",the water can get out ofit.If stepping on the button of blow, the handpiece can blow, and you can adjust the water mist volume of handpiece with spray knob. And you can adjust the working air pressure with air pressure knob aside(normal air pressure is 0.25 MPa). Get off the foot controller to get power off,then put the handpiece back on to the rack.

6.4 Low speed handpiece

Take the end of low-speed handpiece out of the rack,and connect it to the low speed handpiece,the operation steps are the same with those of high speed handpiece. You can adjust the water mist volume of handpiece with spray knob,and you can adjust the working air pressure with air pressure knob aside(normal air pressure is 0.3MPa). Get off the foot controller to get power off, then put the handpiece back on the rack.

Notes

We'd better raise the water or air volume gradually to get the ideal performance.do not idle the handpiece or start with over airpressure.otherwise it may shorten lifespan of the handpiece.

6.5 3-ways syringe

This dental unit is equipped with 2 pieces of 3-ways syringe.When using the 3 way syringe,press the corresponding button to eject water(air), while pressing the song button can be sprayed

6.6 Stronge/Weak saliva ejector

The stronge ejector or weak ejector will start at once when they get away from the rack, on the

other hand,they will stop when they get back to the rack

Air Master Switch

Water switch



Press:On

Pull:Off



Pull:Choose A bottle(A bottle for disinfection water)

Press:Choose B bottle(B bottle for pure water)

Swich under the Touch sense instrument tray

7.Pure water bottle

7.1 One bottle & only for R7 implant model

This bottle is equipped for use for 3-ways syringes and handpieces of instrument tray.turn the air switch onto ON. And turn on water switch,then it will change to get the pure water ready. Close the air supply witch before adding pure water into the water bottle. Turn the switch to OFF side if needing to use running water.

7.2 Two bottles(only for R7 and R7 floor mounted)

Compare with single-bottle dental unit,Double bottle dental unit has two switches for two bottles for two bottles'control,usage is the same as said above,the single one,the switch at the side of A, means the bottle A is ready to suppl water. Turn to B side,means the bottle B is ready to supply water(Front for A bottle,behind for B bottle)

7.3Function of cleaning and disinfection(only for R7 and R7 floor mounted)

This is the function of cleaning and disinfecting handpiece water tube. Operation steps described as below: A) preparation jobs

Add pure water into bottle A,and add bottle B the disinfectant liquid.(advice of 500mg/L disinfectant liquid including chlorine). Shift the master switch to ON side to get the air and water supply ready. Turn the water switch to the side A

B) Flush and disinfection

Take the hadpiece connective end out of the rack and put it over the spitoon,turn on the disinfection switch at the instrument tray,then water spouts out from the end of handpiece connective end.and you can control the water volume from the end with valve knob beside the instrument tray.(20 to 60 seconds for advice)

Turn the pure water bottle switch to side B,the handpiece can spout out infectant liquid at this time.(2 to 10 minutes of sterilisting the handpiece tube for advice), when turning the switch to A side again, it will flush the handpiece tube again to clean away the disinfectant liquid (we advise 20 to 60 seconds)




C) Turn the disinfection switch to OFF side, Put handpiece connective end back to the rack after finishing flushing the tube.


Notice

There is a air pressure reducing valve in the man machine box for controlling its level during 0.2-0.25MPa, but it is set before delivery to user, please do not adjust the valve except the professional specialized man.

8.Usage and operation direction



R7 main panel



R7 assistant panel

Symbol	Meanings	Symbol	Meanings
	Seat rising up		Setting button
	Seat getting down		X-ray film viewer button
	Backrest rising up		Dental Lamp
	Backrest getting down		Heating
	Shock position		Tumbler filling
	Entry/exit position		Flushing spittoon
	Instrument tray rising up		Comprehensive key for Flushing,lighting turn off and chair pose recovery
	Instrument tray getting down		Preseted Chair pose 1
	Election of doctor		Preseted Chair pose 2
	Chair pose for patient's spitting		Preseted Chair pose 3
	Increase		Reduce

8.1 Tumbler filling

- Keep pressing the setting button for 5 seconds,buzzer sound and the signal lamp light up at the same time.meaning the it is ready for setting now.
- Keep pressing the button of water supply to cup,to flow water to the volume needed,the light turn on at this moment ,then stop pressing, the water stop flowing out.
- Press the setting button again, the buzzer sound again, the signal lights turn down at the same time.This means that water volume in the cup is the setting volume for next time or future.

8.2 Time setting for flushing spittoon

- Keep pressing the setting button for 5 seconds,buzzer sound and the signal lamp light up at the same time.Meaning the it is ready for setting now.
- Keep pressing the button of flushing spittoon,the light will show to "1",after spittoon will flush disinfection water,after 3-5 seconds,unstuck the button, the wash pipe will turn to tap water and show light"2". This for disinfection timing setting.

Notice

- A)Put 500mg/L of disinfectant liquid including chlorine and turn on the master air switch before setting.
 B)Disinfection time can be set for 0 to 60 seconds.
 C)You can skip over this setting step if you do not need this function,or turn on the air master switch sown the main box.

C) Setting time of flushing spittoon according to times of pressing buttons,one pressing to a sound

Pressing times	Buzzer sound	LED Flash	Flushing time
1 times	1 sound	1 times	3 minutes
2 times	2 sound	2 times	15 minutes
3 times	3 sound	3 times	30 minutes
4 times	4 sound	4 times	No time limited

D)Press the setting button again,the buzzer sound again,this signal light turns down,that means that the time of flushing spittoon is remembered by machine.

Notice

The time of flushing spittoon can be set according to setp A),C),D).

8.3 Seat rising up:Pressing this button to rese the seat up,stop pressing to stop rising

8.4 Seat getting down:Press this button to get seat down,stop pressing to stop getting down.

8.5 Backrest rising up:Press this button to rise backrest up,stop press to stop rising.

8.6 Backrest getting down:Press this button to get the backrest down,stop pressing to stop getting it dow.

8.7 Instrument tray rising up: Press this button to rise up the instrument tray,stop pressing to stop rising.(The function only for R7+/R7+CART/R7 implant)

8.8 Insturment tray getting down:Press this button to lower the instrument tray,stop pressing to stop lowering.(The function only for R7+/R7+CART/R7 implant)

8.9 Patient chair initial pose: Press this button,the backrest rise and the seat get down to lowest, motor can only during its limited.

8.10 Shock position:Press this button, the backrest will lay down to lowest(-5°) but the seat rise to the top.

8.11 Comprehensive key for flushing,lighting down and Chair pose recovery: Press this button,backrest rises to top,and the seat does as well, and the light turns down, spittoon is being flushing

8.12 Chair position for patient's spitting: Press this button, backrest rise to top, dental lamp turns off meanwhile,and water is flushing the spittoon, press again,return to original treatmeant pose,light up and stop flushing.(this function only for X3 and X3+)

8.13 Setting button :the button for preparing to set parameter of other function,press this button again mean to lock the parameter.

8.14 Heating: this button is for control of heating water,it owns the protection system of over heating.and the the temperature is set before delivery to user.

Notice

Heating button of R7: the heating button of the assistant control panel control the heater in the main machine

box, the heating button on the instrument tray can control the heater in it.

8.15 Water supply to cup: Press this button, to set the water volume according to the water flowing time only after pressing the button of setting.

8.16 Flushing spittoon: Press this button to set the water volume to flushing according to the water flowing time only after pressing the setting button.

8.17 X-ray film viewer: This button is for turning on or off viewer light.

8.18 Dental lamp, This button is for controlling the dental lamp luminanc weak or stronge, on or off power (please refer to 6.2 operation instruction of dental lamp)

8.19 Election of doctor : Press this button can choose the preseted patient chair position of different dentists, there are 3 dentists options A, B, C, and every dentist can preset 3 chair position, 1, 2, 3, 1) please keep pressing the "setting button" for 5 seconds until its signal light up, 2) press A, B or C which refer to dentist code, 3) adjust the patient chair position to you ideal position by pressing the seat and backrest buttons, set it by press P1


4) Same does the other chair position 2 or 3.

5) Press the setting button again to lock your choice.

Set the temperature, time and display pressure

Press on the control panel "+" or "-" key could be display the "air pressure, heater temperature and time of operation" this function on panel.

A) **Setting the heater temperature:** press "+" or "-" key will display the temperature on panel, after press the setting key about 3 seconds, the buzzer sounds "Di", the temperature display will flashing that could press the "+" or "-" key to adjust the temperature. After setting is ok, please press the setting key again means confirm setting. The flash will stop. (Normally please set the temperature 35 to 45 degree)



B) **Setting dentist operation time:** press "+" or "-" key will display the time on panel (only the number), after press the setting key about 3 seconds, the buzzer sounds "Di", the numbers will display as flashing that could press the "+" or "-" key to set the time. Normally for Anesthesia, whitening, curing, dentist need to set the time on phone and waiting, now R7 dental unit could set on panel to remind our dentist, like the alarm clock. please press the setting key again means confirm setting. The flash will stop. Press the "setting" key again, will start to countdown.

Other functions:

A) we can not press the two engines of backrest and seat to work at the same time

B) Different buttons of engine runing parts can lock each other for considering the safe operation

C) the engine will be locked as soon as the handpiec work (the function that handpiece and engine lock each other)

Notice

- A) Recovery pose of patient chair do not need to be set, it runs up to its limited
- B) Heating button work as function of start or finish, its signal light will b flashing when it heating, and it will keep lighting when finish heating. the temerature is set well before delivering this device to user. do not set again.
- C) Button of X-ray film viewer only work for start or finish lighting, do not need to set anything.
- D) All the settings can be stored automaticly against power failure.
- E) There will not be anthing wrong with this device when power failure happens, please go on to work
- F) The buzzer will sound up at once when switching on this dental unit, the power signal lamp light up, the device examine its own system automaticly at this moment.
- G) You can not go on further operation during setting parameter unless the signal lamp stops flashing.

Warning

Please make sure this dental unit can stretch its arms at utmost. thus to avoid knocking with other objects, which may result in incident.

8.20 Usage and operation direction of foot controller:

A) Movement of foot controller

This switch can control backrest and seat of patient chair to move at four direction, back or front, up or down.

B) Water switch of foot controller

Step on this switch to get water supply to cup

C) Lamp switch of foot controller


step on this switch to get the lamp luminace weak or strong, or power on or off.

D) Blowing switch of foot controller

Step on this button to blow away the scraps in the mouth with handpiece

E) Air supply pedal

Step on this pedal to supply air to handpiece



8.21 Handpiece driving air control

Handpiece driving air control knob is below the instrument tray, their symbols are:



Controlling methods: turning clockwise is for reducing air

8.22 Handpiece cooling water control

Handpiece cooling water control knob is located at the right side of instrument tray, their symbols are:



Handpiece cooling water control methods: turning clockwise is for reducing water

8.23 Handpiece spray control

Handpiece spray control knob is located at the right side of instrument tray. its symbol is above.

Troncolling methods: Turning clockwise is for reducing spray.

8.24 Adjustable rack use and operation

A) The instrument tray rack can adjust the angle.

Specific method: Both side have the screw, please make the screw out about 3~5mm, after use the finger to press the screw that will see the rack can adjust the angle. Make the right angle and lock the screw is ok.



Screw position at both side

B)The rack shell remove and install.

Specific method:Remove,two finger press to centre that can remove the rack shell.

Install:The rack shell put it back is ok,will fix it.



Both sides to prssure

8.25 USB charger

Have USB charger to charge the mobile phone on cabinet left side.



8.26 Assistant control drawer

The assistant control have drawer to put the disposable.



8.27 Spittoon lock(only suit for R7 implant and R7 floor mounted)


Spitton can rotate 180°. In parallel to the seat cushion direction as a benchmark, swing chair frame will lock in the direction to the seat cushion.Parallel with the chair frame or to another direction, chair frame is unlocked.



8.28 Electrocardiogram equipment and Intraoral camera(only suit for R7 implant)

Electrocardiogram equipment:The doctor can observe the patient's heart electrical activity

Intraoral camera:The doctor can real-time video inspections need to check the parts.



8.29 Multi-function stainless steel plate

3 stainless steel plate,different size,can remove,clean and disinfection.



7.Maitenance and preservation

For daily maintenance, Please keep this device clean, keep passage of tube system clear,examin that there is no leaking of water or air.the rotatable parts should be applied to lubricant oil on schedule to keep it at well performance

7.1 water filter

There is a water filter at the end of water-in tube connecting in the base box, this is for blocking impurity in water to get into the device, user should clean or exchange this filter on schedule to make its passage clear.

7.1.1 In general,if there is one of the condition listed below, means the time to clean or exchange the water filter.

- A)This dental unit has been used for over a year
- B)Pressure loss from this filter is over 0.1MPa
- C)This filter is indefected
- D)Filtered water turn out muddy

7.1.2 Cleaning or exchanging steps

- A)Shut down water and power supply
- B)Open the front base cover,rotate down in anticlockwise the crystal cup under the filter.
- C)Unscrew the fixed screw in the midle of the filter, tear it off to clean or exchange
- D)Fix back the new or clean one in reverse steps

7.1.3 Right material and size for the filter

Water filter:its filtering degree of accuracy is 40um, it is sintered with bronze particle

Warning

Can not overexert while installing the water filter, avoid crushing the filter. Keep the crystal cup in well sealed.

7.2 Air filter regulator

This dental unit is equipped with an air filter regulator at the air-in end of the air supply tube in the front base box, in order to ensure the air supply is clean, stable and dry, as well as reducing the air pressure meanwhile. It can also filter impurity and water in the air. So there must be some wasted water gathered in the filtering cup, we should clean it after a long time of use.

7.2.1 If there is one of conditions listed as below, we must empty the water in the air filter regulator.

- A) The dental unit has been used for over a week
- B) Water volume in the bottle reaches up to 3/4 of total volume
- C) The water in the bottle turns its color (not transparent water)

7.2.2 Emptying water steps of air filter regulator.

- A) Shut down electricity and air supply
- B) Open the front base cover, release the rest air in the device by pressing 3-way syringe,
- C) Take down the water bottle by turning clockwise, unscrew the air filter regulator to get rid of the water inside

7.2.3 The size and material match with this air filter regulator

Air filter regulator: filtering accuracy is 5µm, material is polyethylene fibre

7.3 Spittoon

Spittoon must be cleaned on schedule, the filter net at the drain must be cleaned every day, to prevent blockage there.

7.4 Patient chair

Clean and sterilize the patient chair in time with disinfectant liquid or neutral detergent.

7.5 Dental Lamp

Please maintain it referring to the user manual of dental lamp

7.6 Filter net of saliva ejector

The two filter nets are inserted at the joint parts between the ejector tube and the main machine box, they should be cleaned on schedule (better for every week), which usually hold some dirty scraps when ejecting water or air pass through them.

8. Malfunction removal

Malfunction Description	Malfunction Causes	Removal methods
1.Powerless handpiece	1.Bad bur	Contact the handpiece seller
	2.Bearing fault	
	3.Low driving pressure	Check the air pressure
2.Powerless saliva ejector	1.Air or water pressure not high enough	Raise the air or water pressure
	2.Joint part not sealed well	Exchange the seal ring
3.Leaking water from the tap to cup or the one for flushing spittoon	3.Magnetic valve gets blocked	Clean the Magnetic valve core
4.Driving air pressure not stable	1.Filter gets blocked	Clean the filter
	2.The tube gets twisted	Stretch the tube
	3.Air supply not open	Open the air supply
5.No water mist gets out from handpiece	1.Water supply gets blocked	Contact handpiece seller
	2.The water control knob get wrong ways	Refer to user manual
6.Handpiece water control fail	Valve for controlling water break down	Refer to user manual
7.Handpiece spouts out water and bubble	1.Shim of handpiece breaks down	Contact handpiece seller
	2.Handpiece tube gets broken	Examine and exchange
8.Handpiece sprays badly	1.Handpiece spray nozzle gets blocked	Contact handpiece seller
	2.The water control knob get wrong ways	Turn control knob to right level
9.Leaking water at the joint part of handpiece	1.Shim of high speed handpiece breaks down	Contact handpiece seller
	2.Too much water in the air filter	Empty the filter every day
10.No water spout out from 3-way syringe	1.Its water tube gets blocked	Block the spray nozzle,and press water and air control button
	2.The water control knob gets the wrong ways	Adjust the control knob again
11.Air mingled with water gets out while 3-way syringe spurts air	1.Too much water in the air filter	Empty water in time, examine to exchange bad parts, or install again
	2.Get something wrong with air valve core	
	3.Wrong installation of spray nozzle	

12. Leaking air when handpiece stops working	Rack valve gets loose	Adjust the rack valve
--	-----------------------	-----------------------

9. Waste treatment

No handling these dental unit waste by ways of household waste done They must be entirely disinfected,sterilized and purified before being dismantled or being thrown away

Part names	Hazardous substances or element					
	Lead (Pb)	Hydrargyrum (Hg)	Cadmium (Cd)	Hexavalent chromium (Cr6+)	Polybrominated biphenyl (PBB)	Polybrominated DiphenylEthers (PBDE)
Power wire with plug	×	○	○	○	○	○
Power switch	×	○	○	○	○	○
Voltage transformer	×	○	○	○	○	○
Fuse	×	○	○	○	○	○
Over temperature Protector	×	○	○	○	○	○
Electric Motor	×	○	○	○	○	○
PCB	×	○	○	○	○	○

○: This mark means the hazardous material content in the homogenous parts is under the standard regulated in SJ/T 11363-2006 , the REQUIEMENTS FOR CONCENTRATION LIMITS FOR CERTAIN HAZARDOUS SUBSTANCES IN ELECTRONIC INFORMATION PRODUCTS

×: This mark means that hazardous substance content in a certain part of this homogenous parts is beyond limits required by SJ/T 11363-2006

Explanation for the parts with lead element

1. Power wire with plug: most parts of power wire with plug or switch including the metal cover of the connector, ends of the pins
2. High temperature solder contain 85% lead, which is adopted to welding the voltage transformer.
3. The resistance in the fuse
4. The over temperature protector resistance contain lead

5.The resistance in the electric motor contain lead

6.PCB bonding pad contain lead

10.Taboos, warnings and notices

1.Read and master the user manual completely before working on this device

2.Operate according to the warnings and instructions of this device.

3.Power wire must the requiemment of standard,the ground wire must be connected to ground well.

4.Never heating water without water supply,otherwise , that must break down heater.

5.Max load of this instrument tray is 2.5KG, pushing, pulling or pressing this try should be soft instead of overexerting.

6.Air pressure of this dental unit has been set down before delivery to user, do not adjust it by yourself unless the careerman.

7.Action limit of this dental unit has been locked,do not adjust it by yourself unless the careerman.

8.Shut down the master switch of water,air,electricity before getting off work.

9.If the foot controller adopts the round one with switch you should turn on this switch first,then step on this pedal to supply water to handpiece.

10.For this dental unit users, the dental clinic or other setups. amalgam separator must be equipped with this dental unit at the end of waste discharge system.

11.Before maintenance, we must empty the water and air first , and then shut down power supply.Open the water,air and electricity supply before using it.

12.Handle and spray nozzle of 3-way syringe should be put in a sterilizer bag first before sterilization, and do high sterilization at 135°C (2 bar), for less than 15 minutes.

13.For sterilization of high speed handpiece and low speed handpiece, please refer to the handpiece user manual, Sterilization time should be over15 minutes.

14.Use wet colth with neutral cleaner to clean the device shell, to ensure the surface smooth and shining,and to keep it away from corrosion, Use soft cloth with soap water to clean PU and ABS parts, Use wet colth with soap water or wax spray to clean the baking finish of metal parts.

15. Mind to avoid damaging the optic fiber part when you use a fiber optic handpiece,keep it away for a distance to avoid touching the waste produced from operating the handpiece,you can use a optic shelter 5 seconds before operating handpiece.

16.Get rid of any dirty on tools at once, take off handpiece and clean with a cloth alcohol after finishing using it.

17.Mind not to let the light from handpiece optic fibre shining on patient' s eyes! which may hurt patients' eyes. like those of cataract. in general, these rays do not damage eyes constantly,but it can result in temporary blindness.

18.Maximum load of patient chair is 135kg.overloading may cause it fail to work.

19.when using extenal equipment to do dental implant operation together with this dental unit, but if no needs of dental unit to supply electricity to extenal equipments, please shut down the power supply to dental unit. to avoid bringing accident due to triggering any switch mindlessly.

20.Take off handpiece burs only after the handpiece stop working completely, otherwise the bars holder may get broken or fall down to hurt someone.

21.Use high quality bars and right size of screw bolt.

22.Check the bars holder everyday before using handpiece if the holder hold bars firmly.

23.Try to pull out the bars of high speed handpiece to check if it is installed well after it is

exchanged.

24. Diameter of bars should be between 1.59 to 1.60mm (regulation by ISO1797 Class 3), Max length is 25mm (regulation by ISO6360-1)

25. Only after being equipped with bars and dental restoration tools, high speed handpieces can be used

26. Do not press the off button of bars while the handpiece is working, the friction force between the button and the impellers of the air driving engine may burn the handpiece head due to overheating.

27. The tissues of the patient's mouth (tongue, cheeks, lips), etc., must be protected by appropriate means (with mirrors, etc.) so as not to touch the buttons.

28. Put on gloves to exchange dental lamp to avoid to be burnt by the lamp, do not touch lamp directly with bare hands.


29. Do not touch the electronic parts on the PCB with bare hands or metal tools

30. Air supply pressure is between 0.55MPa and 0.8MPa, and air supply pressure is 0.2-0.4MPa, the water should be clean enough.

31. Examine air and water tubes system on schedule if there is leaking


32. Mobile or portable high frequency communication tools may disturb medical electronic equipment. It is not allowed to use cell phone or other mobile communication tools in hospital or clinic. If you need to use external high frequency surgical equipment, please do not hang it onto this dental unit.

11. Electrical schematic diagram



12. Water-air layout diagram

R7 air layout diagram



1. Crystal tubes are for air layout system
2. Blue PU tube are for air layout system
3. PVC tubes (25MM of diagram) are for waste outlet