

## JT-25 Rotatable Gas Lights Exploding Diagram(NEW)

# JT-25 Rotatable Gas Lights Exploding Diagram(NEW)

Serial No.	Item No.	Q' ty
1	Working Tube Kit	1
2	Kindling Tube Kit	1
3	Flame adjustment Ring(L)	1
4	Flame adjustment Ring(S)	1
5	Clapboard	1
6	Wave pad	1
7	Rotatable Screw	1
8	O Ring	6
9	Base	1
10	Working Tube Valve	1
11	Backplane	1
12	Rubber Support	4
13	Back Valve	1
14	O Ring	1
15	Front Valve	1
16	O Ring	1
17	Air-supplier Tube	1
18	Kindling valve	1
19	The total valve	1

This is the latest ,energy-saving gas lights with best configuration . Patent expected to get approval soon. It comes with below-mentioned specifications:

- 1, It can rotate 30 degree freely from left or right. It works well both for left-handed and right-handed
- 2,Kindling Tube and Working Tube work separately . They both have their own controlling system .It make this gas light save energy .and with single-tube and double-tube function
- 3,Distinguished controlling system ensures the gas light never leak gas during working

#### Operation Manual

1. Pls turn off 10/18/19 valves before connect with gas
- 2, Connect the gas supplier with Air-supplier tube, turn on gas-supplier valve. After confirming that no gas leak , turn on the air-supplier valve 19. then turn on the kindling valve 18 slowly. Make the fire on ,turn on the working tube valve 10, then adjust the working flame
- 3,Control the flame strong/weak by controlling flame adjusting Ring 3/4
- 4, You can turn off the kindling valve to make it working in single tube
- 5,Remember to turn off fires and all valves after duty off
- 6,This gas lights is suitable for Bottled or pipeline Nature Gas/Coal Gas/Oil-gas

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