

AIXIN MEDICAL

USER MANUAL

EM-BX1

DUST CABINET

Description:

The EM-BX1 Dust Cabinet features wide operating view, strong cold lighting system, wide hand inlets, dust collection system, and the optional wooden arm rests. With this cabinet, dental technicians no longer have to do the irritating desk top cleaning work, and it can get rid of dusts, fragments and gas, which greatly improves working quality and working environment of dental technicians.

Specifications:

Voltage	220V / 110V 50Hz / 60HZ \pm 10%
Power	<8W
Optional part	Dust Collector
Arm Rest	Removable
Dimensions	39 × 34 × 26cm
Package size	44 × 39 × 31cm
Weight	15kg