

What is invisible aligner?

Invisible orthodontic treatment is based on the professional integration of 3D printing technology, CAD/CAM and other sophisticated technology. alinger is made of imported oral special polymer film (PETG) which is pharmaceutical grade raw materials.

Working Principle: After wearing it, it relys on the resilience force which generated by material elastic deformation to straighten teeth. Thereby moving teeth and reach the correction purposes from malocclusion.

How does invisible aligner work?

Invisible orthodontics includes a series of transparent aligners,and each aligner is exactly designed so that the teeth make tiny movement. To wear each stage accordingly makes teeth moved gradually and reached the ideal position.



Professional digital dental provider

Difference between invisible aligners and traditional brakets

Compared with the traditional brakets , invisible aligners can not only achieve good results, but it has incomparable advantages such as invisible, comfortable, beautiful and suitable for most people who need treatment.

The 3 steps of invisible aligners

- 1 The dentists collect patients’ silicone rubber impressions, intraoral pictures, and patients form, then deliver to invisible aligners
- 2 Invisible will send animations and design sheet to dentists within 7 working days.
- 3 After the dentists confirmed animations and paid, will deliver invisible aligners within six working days.



Invisible Aligner Indications

- Mild to moderate crowding
- Diastema
- Dental crossbites
- Deep overbite
- Extraction cases
- Relapse
- Patients with periodontitis
- Orthognathic surgical preoperative and postoperative treatment
- Adjustment of the end of lingual orthodontics



Invisible Aligner

Discover beauty witness your smile mircle

FAQ

Am I suitable for treatment with invisible alingers?

Invisible treatment is suitable for most patients who need orthodontic treatment. If you need to make invisible treatment, you had better consult clinical dentists.

How long invisible orthodontic treatment will take?

Treatment time depends on patient's malocclusion situation. It normally takes 6-18 months.

How does the orthodontic effect by invisible aligners?

A satisfied and effective treatment results will be achieved as long as patients treated by professional clinical dentists and wear invisible aligners carefully according to instructions.

What should invisible orthodontic treatment pay more attention?

The aligners can be removed while eating, brushing teeth or on social occasions. Patients should remain wearing at least 20 hours daily. Careful cooperation is tightly related to good treatment effect.

Advantages

Invisible	Effective	Comfortable
Customized	Predictable	Removable

