

# Examination Instrument Kit

Tray,Mirror,Explorer,Pliers,Bibs



## Matters need attention:

The Stamato-kit Disposable is for one-time dental clinic use,it was sterilized by Epoxy Ethane. The Stamato-kit Disposable is valid for 2 years.It can't be used if the package is broken,things get rust,the kit out of the expiration date.

## Using method:

Open the plastic bag,using one time for one person only.

## Storage:

Please put this product in the cool dry in door without corrosive gas,away from the file.

This product is composed by:Tray,Mirror,Explorer,Pliers,Bibs

---

### Mirror:

Size:L17cm /Mirror Diameter:2.5cm



### Pliers:

Size:L14.8cm



### Explorer:

Size:L17.5cm,D1.3cm



### Tray:

Size:19\*11.5\*2cm



### Bibs:

Size:L44\*D33cm

